

# Titan Taq DNA Polymerase

Cat. No.	Pack Size	Conc.
BTS0001	50 U SAMPLE	5 U/μl
BT10101	500 U	5 U/μl
BA00102	1000 U	5 U/μl
BA00103	2000 U	5 U/μl

## Storage & Shipping:

Store at -20°C, shipping at room temperature.

Shipping and temporary storage for up to 1 month at room temperature has no detrimental effects on the quality of Titan Taq DNA Polymerase.

## Reagents Provided:

- **Titan Taq DNA Polymerase**
- **10x Reaction Buffer 1 (Mg<sup>2+</sup> free):** 800 mM Tris-HCl, 200 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 0.2% w/v Tween-20.
- **25 mM MgCl<sub>2</sub>.**

## Description:

Titan Taq is a highly processive, thermostable DNA polymerase. Due to its genetic modifications Titan Taq has an enhanced stability at room temperature with no activity loss for up to 1 month.

The enzyme has 5'→3' polymerization-dependent exonuclease replacement activity but lacks 3'→5' exonuclease activity.

## Source:

Purified from an *E. coli* strain that carries an overproducing plasmid containing a modified gene of *Thermus aquaticus* DNA Polymerase.

## Storage and Dilution buffer:

50% glycerol (v/v), 20 mM Tris-HCl pH 8.7 at 25°C, 100 mM KCl, 0.1 mM EDTA and stabilizers.

## Quality data:

The enzyme is free of nicking and priming activities, exonucleases and non-specific endonucleases. SDS/PAGE - 95 kD band, >98% pure. Activity and stability tested via thermo-cycling. The error rate per nucleotide per cycle is ~ **2.5 x 10<sup>-5</sup>**; the accuracy is ~ 4 x 10<sup>4</sup>. Estimated half-life at 95°C is 1.5 hours.

## Unit definition:

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nmol of dNTPs into an acid-insoluble form in 30 minutes at 74°C.

#### Recommended PCR reaction mix:

Component	Volume	Final conc.
Titan Taq (5 U/μl)	0.4-1.0 μl	0.02-0.05 U/μl (2-5 U)
10x Buffer 1	10 μl	1x
25 mM MgCl <sub>2</sub>	6-10 μl	1.5-2.5 mM
20 mM dNTP mix	1 μl	200 μM
Primer Forward (10 pmol/μl)	1-3 μl	0.1-0.3 μM
Primer Reverse (10 pmol/μl)	1-3 μl	0.1-0.3 μM
DNA template	5-20 μl	5-100 ng/μl
H <sub>2</sub> O PCR grade	Up to 100 μl	
<b>Total</b>	<b>100 μl</b>	

#### Recommended PCR cycles:

Cycle step	Temp.	Time	Cycles
<b>Initial denaturation</b>	<b>95°C</b>	<b>3-5 min</b>	1
Denaturation	95°C	30-60 s	26-35
Annealing	50-68°C	30-60 s	
Elongation	72°C	1-4 min	
Final elongation	72°C	5-10 min	1

**IMPORTANT:** Annealing temperature should be 2-6°C lower than the primer melting temperature. Elongation time should be ~1 min/1 kb.

#### Safety warnings and precautions:

This product is designed for research purposes and *in vitro* use only. According to common laboratory safety practice, it is recommended to wear protective clothing, gloves and safety glasses. Please refer to [www.bioatlas.com](http://www.bioatlas.com) for Material Safety Data Sheet of the product.

*Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required.*