

Atlas 10x Taq Buffer

Cat. No.	Pack Size
BA00151	1.8 ml

Storage & Shipping:

Store at -20°C, shipping at room temperature.

Description:

Atlas 10x Taq Buffer is a ready-to-use buffer to be used for Atlas Taq, Atlas HotTaq, Atlas RedTaq and Atlas LongTaq enzymes in PCR reactions. Does not contain MgCl₂.

Composition:

100 mM Tris-HCl (pH 8.8 at 25°C)

500 mM KCl

0.8% Nonident P40

Safety warnings and precautions:

This product is designed for research purposes and *in vitro* use only.

According to common laboratory safety practice, it is recommended to wear protective clothing, gloves and safety glasses. Please refer to www.bioatlas.com for Material Safety Data Sheet of the product.

Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required.