

Atlas 25 mM MgCl₂

Cat. No.	Pack Size
BA00153	1.8 ml
BA00154	5 x 1.8 ml
BAS0022	1.8 ml SAMPLE

Storage & Shipping:

Store at -20°C, shipping at room temperature.

Description:

Atlas 25 mM Magnesium Chloride (MgCl₂) can be used for optimization of magnesium ion concentration in PCR reaction.

Safety warnings and precautions:

This product is designed for research purposes and *in vitro* use only. According to common laboratory safety practice, it is recommended to wear protective clothing, gloves and safety glasses. Please refer to www.bioatlas.com for Material Safety Data Sheet of the product. *Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required.*