

Atlas ClearSight Agarose Tablets

Safety Data Sheet

Bioatlas Ltd Posted
info@bioatlas.com

Date: 04.10.2018
Print Date: 06.10.2018
Version 1
Revision Date 13.08.2019

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifiers

Product name:	Atlas ClearSight Agarose Tablets
Product Number:	BA50202
Brand:	Bioatlas
EC Substance Name:	Not applicable
EC number:	Not applicable
CAS reference number:	Not applicable
REACH reference number:	A registration number is not available for this substance as the substance or its uses are exempted from registration, the annual tonnage does not require a registration or the registration is envisaged for a later registration deadline.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses: For laboratory and manufacturing use only.

1.3 Details of the supplier of the safety data sheet

Company: Bioatlas
Riia 181A
Tartu, 51014
Estonia
Phone: +372 5562 4867
e-mail: info@bioatlas.com

1.4 Emergency telephone number

Emergency Phone: +372 5562 4867
* Within business hours: Monday to Friday, 8am to 5pm

SECTION 2: Substance or mixture possible hazards classification

2.1 Classification according to EC Directives 67/548 / EEC (DSD) and 1999/45 / EC (DPD):

Not a hazardous substance or mixture according to EC Directives 67/548 / EEC (DSD) and 1999/45 / EC (DPD)

2.2 Classification according to Regulation (EC) No 1272/2008 (CLP):

Signal word: Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008

H-Statements: Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008

Precautionary statement(s)

P280	Wear protective clothing (gloves, eye or face protection)
P301+ P312	IF SWALLOWED: call poison center or doctor if you feel un well
P302+ P352	IF ON SKIN: Wash with plenty of water
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308 + P313	IF exposed or concerned: Get medical Supplemental advice/ attention.

Hazard Statements None

2.3 Other hazards:

None

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

2.4 Additional hazards for humans and the environment:

2.4.1 Possible hazardous physicochemical effects:	Unknown
2.4.2 Possible adverse effects on humans and possible symptoms:	No data available
2.4.3 Possible harmful effects:	No data available

Restricted to professional users.

SECTION 3: Composition / information on ingredients:

3.1 Chemical characterization:

Powder solution with inorganic and organic components.

Product name: Atlas ClearSight Agarose Tablets

3.2 Components

Name	CAS No.	EC No.	Proportion %	Classification according to 67/548/EEC and 1999/45 / EC	Classification according to regulation (EC)No1272/2008
Atlas ClearSight DNA Stain *	-	-	<1%	-	-
Agarose	9012-36-6	232-731-8	>99%	-	-

*No components need to be disclosed according to the applicable regulations.

SECTION 4: First aid measures

4.1 General information:



If skin irritation or other symptoms persist, consult a physician. Show this safety data sheet to the doctor in attendance. Remove contaminated clothing and shoes, carefully clean before reuse.

4.2 After inhalation:

Move to fresh air. If not breathing, give artificial respiration. Get medical attention if symptoms occur.

4.3 After skin contact:

Wash off immediately with plenty of water for at least 15 minutes. If skin irritation persists, call a physician.

4.4 After eye contact:



Eye contact: Remove contact lenses. Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Obtain medical attention.

4.5 After swallowing:

Never give anything to an unconscious person. Do not induce vomiting. Rinse the mouth, spit out the fluid and drink plenty of water. Get medical attention if symptoms occur.

4.6 Notes to physician:

Treat symptomatically.

4.7 Most important symptoms and effects, both acute and delayed.

No data available.

SECTION 5: Firefighting measures

5.1 Suitable extinguishing media:



This product is not flammable. Use extinguishing agent suitable for type of surrounding fire: water spray, dry chemical or carbon dioxide, fire extinguishing foam, extinguishing powder.

5.2 Extinguishing media which must not be used for safety reasons:

No data available.

5.3 Special hazards caused by the product, combustion products or resulting gases:



Thermal decomposition can lead to release of irritating gases and vapors.
Hazardous Combustion Products: Carbon oxides.

5.4 Special protective equipment for fire-fighting:



Wear protective clothing. Wear self-contained breathing apparatus for firefighting if necessary.

5.5 Additional information:

No data available.

SECTION 6: Accidental release measures

6.1 Personal precautions:

Avoid contact with skin, eyes and clothing, wear appropriate protective equipment. Avoid dust formation. Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust. For personal protection see section 8.

6.2 Environmental precautions:

Should not be released into the environment.

6.3 Methods for cleaning up:

Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

6.4 Reference to other sections.

For disposal see section 13.

SECTION 7: Handling and storage

7.1 Handling:

7.1.1 Information for safe handling:

Wear personal protective equipment. Avoid contact with skin and eyes. Avoid ingestion and inhalation. Avoid formation of dust and aerosols. Avoid exposure - obtain special instructions before use.

Provide appropriate exhaust ventilation at places where dust might form. For precautions see section 2.2.

7.1.2 Information on fire and explosion protection:

See section 5.

7.1.3 Handling rules:

Handle in accordance with good industrial hygiene and safety practice.

7.2 Storage:

7.2.1 Technical measures and storage conditions:

Always keep in containers of the same material as the original.

7.2.2 Packaging materials:

Packing materials must be tested for durability before use.

7.2.3 Requirements for storage rooms and containers:

Keep container tightly closed in a dry and well-ventilated place. If possible, keep in the original container. Do not use food containers because of the risk of confusion. Containers which are opened must be carefully resealed.

7.2.4 Further information on storage conditions:

Recommended storage temperature 4-30°C. Protected from light.

7.3 Specific end use(s):

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Components with workplace control parameters:

Contains no relevant quantities of substances with occupational exposure limit values that have to be monitored at the workplace according to TRGS 900.

8.2 Exposure controls

Appropriate engineering controls:

Handle in accordance with Good Laboratory Practice. Wash hands before breaks and at the end of workday.

8.2.1 Personal protective equipment



Respiratory protection:

In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards. Half masks according to EN 140 or full masks according to EN 136, with filters according to EN 143-P1.



Skin protection:

Handle with suitable protective gloves. Gloves must be inspected prior to use. Damaged and worn protective gloves should be replaced immediately. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use

in accordance with applicable laws and Good Laboratory Practices. Wash and dry hands. The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it. Completely unsuitable are gloves made of fabric or leather. The following gloves are suitable:

For full contact:

Material:	Minimum layer thickness	Breakthrough time:
Nitrile rubber	0.35 mm	≥ 480 min
Latex	0.35 mm	≥ 480 min

For splash contact:

Material:	Minimum layer thickness	Breakthrough time:
Natural rubber	0.5 mm	≥ 120 min
Latex	0.5 mm	≥ 120 min
Polyvinyl chloride	0.5 mm	≥ 120 min
Nitrile rubber	0.2 mm	≥ 30 min



Eye protection:

Safety glasses with side-shields conforming to EN166 U



Body protection:

Long sleeved clothing. Special body protection generally not required.

General protective and hygienic measures

Do not eat, drink or smoke during working hours. Keep away from food and drink. Avoid contact with eyes and skin. Remove contaminated and soaked clothing immediately. Wash hands before breaks and after work.

8.2.2 Limitation of environmental exposure:

Prevent further leakage or spillage if safe to do so.

8.2.3 Limitation of end-user exposure:

Avoid formation of dust and aerosols. Avoid breathing vapors, mist or gas. Remove all sources of ignition.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

- | | |
|---|-------------------------------|
| a) Appearance Form: | White, powder in tablet form. |
| b) Odour: | No data available |
| c) Odour Threshold: | No data available |
| d) pH: | No data available |
| e) Melting point/freezing point: | No data available |
| f) Initial boiling point and boiling range: | No data available |
| g) Flash point: | No data available |
| h) Evaporation rate: | No data available |
| i) Flammability (solid, gas): | No data available |
| j) Upper/lower flammability
or explosive limits: | No data available |
| k) Vapour pressure: | No data available |
| l) Vapour density: | No data available |
| m) Relative density: | No data available |
| n) Water solubility: | Soluble in water |
| o) Partition coefficient: | No data available |
| p) Auto-ignition temperature: | No data available |
| q) Decomposition temperature: | No data available |

r) Viscosity: No data available
 s) Explosive properties: No data available
 t) Oxidizing properties: No data available

9.2 Other safety information: No data available

SECTION 10: Stability and reactivity

10.1 Reactivity

No data available

10.2 Chemical stability

Stable under recommended storage conditions.

10.3 Possibility of hazardous reactions

None under normal processing.

10.4 Conditions to avoid

High pressure and temperatures. Stable under regular laboratory conditions. See section 7

10.5 Materials to avoid

Strong acids, oxidation and reducing agents.

10.6 Incompatible materials

Strong oxidizing agents

10.7 Hazardous decomposition products

Other decomposition products

No data available.

In the event of fire: see section 5

SECTION 11: Toxicological information

11.1 Information on toxicological effects

11.1.1 The acute toxicity test

Component	Parameter	Value	Species
Atlas ClearSight DNA stain	LD ₅₀ oral	1000-10000 mg/kg	Kunming mice

Other components: No data available

11.2.2 Corrosive and irritant effects

Skin corrosion/irritation

Irritant effect

Serious eye damage/eye irritation

Irritant effect

11.2.3 Respiratory or skin sensitization

No data available

11.2.4 Aspiration hazard

No data available

11.2.5 Carcinogenicity, mutagenicity and toxicity to reproduction mutagenicity:

Component	Parameter	Value	Cell culture/ Species	Method	Remarks
<i>Atlas ClearSight DNA stain</i>	<i>In vitro</i> bacterial mutation test	0,5-5 mg/plate	Salmonella typhimurium TA97/98/100/102	AMES-test	Non-mutagenic
<i>Atlas ClearSight DNA stain</i>	<i>In vivo</i> micronucleus test	1000-5000 mg/kg	Kunming mice	The mouse bone marrow micronucleus test	Negative
<i>Atlas ClearSight DNA stain</i>	<i>In vitro</i> mammalian chromosomal aberration test	31,2-5000 µg/ml	Chinese hamster ovary line.	<i>in vitro</i> mammalian cell chromosome aberration detection system	Negative

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

Reproductive toxicity:

No toxicological data are available on the components of the product

11.2.6 Additional toxicological information:

No data available

11.3 Practical experience

No data available

11.4 General remarks:

When used and handled according to instructions, the product does not cause harmful effects according to our experience and current information.

SECTION 12: Ecological information

- 12.1 Toxicity** No data available
- 12.2 Persistence and degradability** No data available
No data on biological and abiotic degradation are available for the components of the product
- 12.3 Bioaccumulative potential** No data available
- 12.4 Mobility in soil** No data available
- 12.5 Results of PBT and vPvB assessment**
This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher
- 12.6 Other adverse effects** At working concentration does not pose a risk to aquatic life
- 12.7 General remarks:**
When used and handled according to instructions, the product does not cause harmful effects according to our experience and current information.

SECTION 13: Disposal considerations

- 13.1 Product:**
The allocation of a waste code number according to the European waste catalog should be done in consultation with the regional waste disposal company.
- 13.2 Packaging:**
Residues in packages should be removed, preferably by rinsing with water, and after complete emptying in accordance with the regulations for waste disposal. Packaging which is not completely emptied must be disposed of in the form as determined by the regional waste disposal company.

SECTION 14: Transport information

- 14.1 UN number** Not applicable
- 14.2 UN proper shipping name** Not applicable
In accordance with GGVSEB, ADR / RID, IMDG / GGVSee, ICAO / IATA-DGR
- ADR/RID, GGVSEB Not dangerous good during shipping
IMDG / GGVSee Not dangerous good during shipping
ICAO / IATA-DGR Not dangerous good during shipping
- 14.3 Transport hazard class(es)**
ADR/RID, IMDG, IATA None

14.4 Packaging group

ADR/RID, IMDG, IATA None

14.5 Environmental hazards

ADR/RID, IMDG, IATA None

14.6 Special precautions for user

No data available

SECTION 15: Regulatory information

15.1 EU regulations:

15.1.1 Chemical safety assessment according to Regulation (EC) No 1907/2006: Substance safety assessments (CSA) according to Article 14 (1) of Regulation (EC) No 1907/2006 (REACH) for the components of the product are neither banned nor restricted

15.1.2 Classification and labeling according to directives 67/548 / EEC (DSD) and 1999/45 / EC (DPD):

Not a hazardous substance or mixture, see section 2

15.1.3 Classification and labeling according to Regulation (EC) No 1272/2008 (CLP):

Hazardous substance or mixture, see section 2

15.1.4 Hazardous components which must be listed on the label.

15.1.5 Pictograms:

None

Signal word: Not hazardous to Regulation (EC) No. 1272/2008

H-Statements: Not hazardous substance or mixture according to Regulation (EC) No. 1272/2008

15.1.6 Special labeling of certain preparations: Neither banned nor restricted

15.1.7 Approvals and / or restrictions of use: Neither banned nor restricted.

15.1.8 Information on Directive 1999/13 / EC (VOC Directive) on the limitation of VOC emissions: Neither banned nor restricted

SECTION 16: Other information

16.1 Training instructions:

Accordance with applicable laws and Good Laboratory Practices.

16.3 Recommended use (s) of application:

For research use only.

16.4 Further information and contact for technical information:

Department issuing data sheet:

Company: Bioatlas
Riia 181A
Tartu, 51014
Estonia
Phone: +372 5562 4867
e-mail: info@bioatlas.com

Bioatlas Ltd

Riia 181A Tartu, 51014 Estonia
Phone: +372 5562 4867
Contact person: Liane Viru
e-mail: liane.viru@bioatlast.com

16.5 Data sources for the preparation of the safety data sheet:

European Chemical Substances Information System (ESIS),

Internet: <http://ecb.jrc.ec.europa.eu/esis>.

Hazardous substance information system of the trade associations (GESTIS), Internet:
<http://www.hvbg.de/d/bia/gestis/stoffdb/index.html>.

Hazardous Substances Database (HSDB) - U.S. National Library of Medicine, Internet:
<http://toxnet.nlm.nih.gov>

Hommel interactive 4.0 - Handbook of Dangerous Goods,

Internet: <http://www.springer.com/dal/home/chemistry>.

CRC Handbook of Chemistry and Physics, 88th Edition, 2007-2008,

Internet: <http://www.hbcnetbase.com>.

16.6 Amended information and reasons for amendment:

Version number:

Version1: 04.10.2018

The reason for the change: Voluntary preparation on the basis of Annex II to Regulation (EC) No 1907/2006 (REACH) and Annex I to Regulation No 1272/2008 (CLP).

The safety data sheet is based on Article 32 (REACH) and is, therefore, voluntarily adopted in accordance with Annex II (EU No. 453/2010) Regulation (EC) No 1907/2006 (REACH).

16.7 Notes:

This information is only intended to describe the safety requirements of the product and is based on the present state of our knowledge. They do not constitute a guarantee for the characteristics of the product described in the sense of the statutory warranty regulations. Please refer to the respective product data sheets for the delivery properties. If the product mentioned in this Material Safety Data Sheet is mixed, mixed or processed with other materials, the data in this Material Safety Data Sheet may not be transferred to the new material, unless otherwise specified

End of Safety Data Sheet