

## Atlas Ethidium Bromide Solution

Cat. No.	Pack Size
BA01901	10 x 1 ml
BA01902	5 x 1 ml
BAS0023	50 µl SAMPLE

### Storage & Shipping:

Storage +4°C protected from light, shipping at room temperature.

### Description:

Atlas Ethidium Bromide Solution (1%) is a concentrated fluorescent staining solution for visualizing ribonucleic acids in agarose and polyacrylamide gels. Ethidium bromide is excited with a standard 302 nm UV transilluminator and emits a red-orange signal.

### Safety warnings and precautions:

This product is designed for research purposes and *in vitro* use only.

According to common laboratory safety practice, it is recommended to wear protective clothing, gloves and safety glasses. Please refer to [www.bioatlas.com](http://www.bioatlas.com) for Material Safety Data Sheet of the product.

*Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required.*