

## Atlas dNTP Set (100 mM each)

Cat. No.	Pack Size	Conc.
BHS0004	4 x 25 µl SAMPLE	100 mM
BH40203	4 x 0.4 ml	100 mM

### Storage & Shipping:

Store at -20°C, shipping at room temperature.

### Applications:

- Standard PCR
- qPCR
- cDNA synthesis
- DNA sequencing
- DNA labeling

### Description:

Separate vials of four dNTPs

Concentrated neutralized solutions

Atlas dNTP Set contains four separate vials of dATP, dCTP, dGTP and dTTP, at a concentration of 100 mM each, in highly purified (>99%) deionized water at pH 7.5.

Atlas dNTP Set is free of RNase and DNase, and suitable for any molecular biology application that requires pure deoxynucleotides, such as PCR, real-time PCR, DNA sequencing, cDNA synthesis and nick translation. Atlas Nucleotide Mixes has been tested for PCR and is compatible for all enzymes provided by Bioatlas.

### Safety warnings and precautions:

This product is designed for research purposes and *in vitro* use only. According to common laboratory safety practice, it is recommended to wear protective clothing, gloves and safety glasses. Please refer to [www.bioatlas.com](http://www.bioatlas.com) for Material Safety Data Sheet of the product.

*Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required.*