

Atlas 10x Taq Buffer with (NH₄)₂SO₄

Cat. No.	Pack Size
BA00152	1.8 ml

Storage & Shipping:

Store at -20°C, shipping at room temperature.

Description:

Atlas 10x Taq Buffer with (NH₄)₂SO₄ is a ready-to-use buffer to be used for Atlas Taq, Atlas HotTaq, Atlas RedTaq and Atlas LongTaq enzymes in PCR reactions. Does not contain MgCl₂.

Composition:

750 mM Tris-HCl (pH 8.8 at 25°C)
200 mM (NH₄)₂SO₄
0.1% Tween 20.

Safety warnings and precautions:

This product is designed for research purposes and *in vitro* use only.
According to common laboratory safety practice, it is recommended to wear protective clothing, gloves and safety glasses. Please refer to www.bioatlas.com for Material Safety Data Sheet of the product.
Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required.