

Atlas PCR water

Cat. No.	Pack Size
BA01201	10 x 1.8 ml
BA01202	25 x 1.8 ml
BA01203	5 x 10 ml

Storage & Shipping:

Store at +4°C, shipping at room temperature.

Description:

Atlas PCR Water is nuclease and DNA contamination free. It is sterile and autoclaved. Manufactured and tested for use as a component in PCR reactions. Validated as well for use with highly sensitive PCR reaction and efficient down-stream cloning applications.

Safety warnings and precautions:

This product is designed for research purposes and *in vitro* use only. According to common laboratory safety practice, it is recommended to wear protective clothing, gloves and safety glasses. Please refer to www.bioatlas.com for Material Safety Data Sheet of the product. *Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required.*